



PRODUCT SPECIFICATION

catalog number:	A21294
product description:	TSH Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human recombinant TSH
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M BBS, pH 8.0
preservative:	None
isotype:	IgG ₁ , Kappa
concentration:	[lot specific] mg/mL
specificity:	Recognizes native and recombinant Human TSH
pairing:	This antibody can be used as the capture antibody in ELISA.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.