



PRODUCT SPECIFICATION

CATALOG NUMBER: A22010

PRODUCT:	Thyroxine (T4) Monoclonal Antibody
IMMUNOGEN:	T4-BSA
FORM:	Protein A column chromatography purified immunoglobulin fraction in 0.015 M KPO ₄ buffer containing 0.85% NaCl and 0.1% NaN ₃ , pH = 7.40
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.