



PRODUCT SPECIFICATION

catalog number:	A22055
product description:	T4 Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1M PBS, pH 7.4
preservative:	0.5 g/L NaN ₃
isotype:	IgG ₁ ,
concentration:	[lot specific] mg/mL
specificity:	Recognizes Natural Human T4
pairing:	Capture Antibody.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.