



## PRODUCT SPECIFICATION

catalog number:	<b>A22511</b>
product description:	Dehydroepiandrosterone-Sulfate (DHEA-S) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Colorless liquid
purification method:	Protein A chromatography; 0.2µm filtered
purity:	>99% by SDS-PAGE
buffer:	20mM sodium phosphate, 0.15M NaCl, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL (based on A <sub>280nm</sub> )
manufacture date:	[lot specific]
storage:	≤-20°C
hazard/biohazard:	There is no known biohazard associated with this product.