



PRODUCT SPECIFICATION

CATALOG NUMBER: A23020

PRODUCT:	Triiodothyronine (T3) Monoclonal Antibody
IMMUNOGEN:	T3
FORM:	Protein A purified IgG in 10mM Phosphate, 150 mM NaCl, pH 7.4
PRESERVATIVE:	0.1% Sodium Azide
AFFINITY CONSTANT:	5.2×10^9
SUBCLASS:	IgG1,kappa
CROSS REACTIVITY:	T3 100% T4 <0.28%
REFERENCE ASSAY:	RIA
STORAGE CONDITIONS:	4°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.