



PRODUCT SPECIFICATION

CATALOG NUMBER: A23070

PRODUCT:	Triiodothyronine (T3) Monoclonal Antibody
IMMUNOGEN:	T3-BSA
FORM:	Chromatographically purified Ig from ascites in 0.05 M PBS, pH 7.4, 0.1% NaN ₃
AFFINITY CONSTANT:	5 x 10 ⁸
SUBCLASS:	IgG
CROSS REACTIVITY:	T4 <1% T2 <1% rT3 <1%
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.