



PRODUCT SPECIFICATION

catalog number:	A23490
product description:	Thyroglobulin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human thyroglobulin
source:	Murine ascites
form:	Liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN ₃
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
specificity:	This monoclonal antibody is specific for human thyroglobulin.
pairing:	This antibody can be used as the capture antibody in a sandwich EIA and pairs with A23480 as the detection antibody.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.