



PRODUCT SPECIFICATION

CATALOG NUMBER: A24010

PRODUCT:	Theophylline Monoclonal Antibody
IMMUNOGEN:	Theophylline-BSA
FORM:	Protein A purified in 10mM Phosphate, 150 mM Sodium Chloride, pH 7.4, with 0.1% Sodium Azide
PURITY:	>90% by SDS-PAGE
AFFINITY CONSTANT:	1×10^8
SUBCLASS:	IgG _{1K}
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.