



PRODUCT SPECIFICATION

catalog number:	A24101
product description:	Tricyclic Antidepressants (TCA) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	TCA-BSA conjugate
source:	Murine
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
reconstitution:	Reconstitute with a neutral pH buffer of choice avoiding extreme high or low pH.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.