



PRODUCT SPECIFICATION

catalog number:	A24140								
product description:	TCA Monoclonal Antibody, Purified								
lot number:	[lot specific]								
immunogen:	TCA-BSA								
source:	Murine ascites								
form:	Liquid								
purification method:	DEAE chromatography								
purity:	>95% by SDS-PAGE and HPLC								
buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, pH 7.2								
preservative:	0.05% NaN ₃								
isotype:	IgG ₁ , kappa								
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)								
specificity:	This antibody is specific for tricyclic antidepressants (TCA)								
cross-reactivity:	<table><tr><td>▪ Amitriptyline</td><td>100%</td><td>▪ Imipramine</td><td>1300%</td></tr><tr><td>▪ Desipramine</td><td>100%</td><td>▪ Nortriptyline</td><td>90%</td></tr></table>	▪ Amitriptyline	100%	▪ Imipramine	1300%	▪ Desipramine	100%	▪ Nortriptyline	90%
▪ Amitriptyline	100%	▪ Imipramine	1300%						
▪ Desipramine	100%	▪ Nortriptyline	90%						
storage:	2-8°C short term (2 years); -20°C long term avoiding multiple freeze thaw cycles.								
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.								