



## PRODUCT SPECIFICATION

catalog number:	A25049																
product description:	Progesterone Monoclonal Antibody, Purified																
lot number:	[lot specific]																
immunogen:	Progesterone-11-BSA																
source:	Murine																
purification method:	Protein A																
purity:	>98%																
buffer:	50mM NaPo4, pH 7.2±0.2, 0.02% Proclin																
isotype:	IgG																
concentration:	[lot specific] mg/mL																
cross-reactivity:	<table><tr><td>▪ 11-a-OH-4-Progesterone</td><td>ND</td><td>▪ Pregnenolone</td><td>&lt;0.01%</td></tr><tr><td>▪ 20-a--OH-Progesterone</td><td>ND</td><td>▪ 5-a-OH-Progesterone</td><td>ND</td></tr><tr><td>▪ Dihydroandrostene</td><td>ND</td><td>▪ 5a-Androstendione</td><td>ND</td></tr><tr><td>▪ dihydrotestosterone</td><td>ND</td><td></td><td></td></tr></table>	▪ 11-a-OH-4-Progesterone	ND	▪ Pregnenolone	<0.01%	▪ 20-a--OH-Progesterone	ND	▪ 5-a-OH-Progesterone	ND	▪ Dihydroandrostene	ND	▪ 5a-Androstendione	ND	▪ dihydrotestosterone	ND		
▪ 11-a-OH-4-Progesterone	ND	▪ Pregnenolone	<0.01%														
▪ 20-a--OH-Progesterone	ND	▪ 5-a-OH-Progesterone	ND														
▪ Dihydroandrostene	ND	▪ 5a-Androstendione	ND														
▪ dihydrotestosterone	ND																
expiration date:	[lot specific]																
storage:	4-8°C																
hazard/biohazard:	There is no known biohazard associated with this product.																