



PRODUCT SPECIFICATION

CATALOG NUMBER: A26050

ANTIBODY TO:	Carcinoembryonic Antigen (CEA)
IMMUNOGEN:	Human CEA
FORM:	Liquid
PURIFICATION:	Purified by A agarose
BUFFER:	10 mM phosphate, 0.15 M sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% sodium azide
SUBCLASS:	IgG _{2a}
STORAGE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.