



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>A26070</b>   |
| product description: | CEA Monoclonal Antibody, Purified                         |
| lot number:          | [lot specific]  |
| immunogen:           | Native Human CEA  |
| source:              | Murine  |
| form:                | Liquid  |
| purification method: | Protein A affinity chromatography                         |
| buffer:              | 0.1 M BBS, pH 8.0   |
| concentration:       | [lot specific] mg/mL                                      |
| specificity:         | Specific for CEA  |
| expiration date:     | [lot specific]  |
| storage:             | 2-8°C   |
| hazard/biohazard:    | There is no known biohazard associated with this product. |