



PRODUCT SPECIFICATION

catalog number:	A26071
product description:	CEA Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Native Human CEA
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	0.1 M BBS, pH 8.0
concentration:	[lot specific] mg/mL
specificity:	Specific for CEA
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	There is no known biohazard associated with this product.