



PRODUCT SPECIFICATION

CATALOG NUMBER: A26150

PRODUCT:	Carcinoembryonic Antigen (CEA) Monoclonal Antibody
IMMUNOGEN:	Human CEA
PURIFICATION:	DEAE column chromatography purified immunoglobulin fraction
BUFFER:	0.015 M Potassium phosphate, 0.85% NaCl, pH 7.40, 0.1% NaN ₃
AFFINITY CONSTANT:	1 x 10 ¹¹ l/mol
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C DO NOT FREEZE