



## PRODUCT SPECIFICATION

catalog number:	<b>A26300</b>
product description:	CEA Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human CEA
source:	Murine ascites
form:	Liquid; filtered in 0.2 µm solution
buffer:	10mM sodium phosphate, 0.80% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This monoclonal recognizes human CEA.
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A26310 (as the conjugate antibody).
notes/applications:	This monoclonal is known to work in lateral flow, TIA and ELISA.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.