



## PRODUCT SPECIFICATION

catalog number:	<b>A27042</b>
product description:	Human CKMM Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human CKMB
source:	Murine
form:	Frozen liquid
purification method:	Protein G chromatography
buffer:	PBS, pH 7.3
preservative:	none
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL
specificity:	Human creatine kinase, MM subunit specific
expiration date:	[lot specific]
storage:	-20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known biohazard associated with this product.