



PRODUCT SPECIFICATION

CATALOG NUMBER: A27060

PRODUCT:	CKMB Monoclonal Antibody
SOURCE:	Murine Hybridoma
IMMUNOGEN:	Purified Human CKMB
PURIFICATION:	Protein G
FORM:	Liquid in PBS, pH 7.3, 0.05% NaN ₃
SUBCLASS:	IgG _{1/K}
STORAGE CONDITIONS:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.