



PRODUCT SPECIFICATION

CATALOG NUMBER: A27370

PRODUCT:	Myoglobin Monoclonal Antibody
IMMUNOGEN:	Human Cardiac Myoglobin
PURITY:	95% by SDS PAGE (matches reference lot) 98.9% by SEC. Distribution of protein bands are comparable to reference lot 149-158- 031804
BUFFER:	20 mM Potassium Phosphate, pH 7.0
PRESERVATIVE:	None
FILTER:	0.2 micro filter
ISOTYPE:	IgG ₁ , kappa
STORAGE CONDITIONS:	-60°C
STATEMENT OF ANIMAL HEALTH:	This antibody is derived from animals with no history of infectious diseases.