



PRODUCT SPECIFICATION

CATALOG NUMBER: A27410

PRODUCT:	Myoglobin Monoclonal Antibody
IMMUNOGEN:	Human Myoglobin
SOURCE:	Immunoglobulin from murine hybridoma
PURIFICATION METHOD:	Protein G purification
BUFFER:	PBS, pH 7.3
SUBCLASS:	IgG _{1/K}
SPECIFICITY:	Myoglobin (17 kD)
STORAGE CONDITIONS:	2-8°C short term, -20°C long term
APPLICATIONS:	ELISA, Western blot, immunoaffinity purification. Pairs with A27420.