



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A27420**

PRODUCT:	Myoglobin Monoclonal Antibody
IMMUNOGEN:	Human Myoglobin
SOURCE:	Immunoglobulin from murine hybridoma
FORM:	Protein G purified immunoglobulin in PBS, pH 7.3
SUBCLASS:	IgG <sub>1/K</sub>
SPECIFICITY:	Myoglobin (17 kDa)
STORAGE CONDITIONS:	-20°C
APPLICATIONS:	ELISA, Western blot, immunoaffinity purification. Pairs with A27410359P.