



PRODUCT SPECIFICATION

catalog number:	A27520
product description:	Myoglobin Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human myoglobin protein
source:	Murine ascites
form:	Liquid
purity:	>90%
buffer:	10mM PBS, pH 7.4
preservative:	0.1% P300
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.