



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A28020**

PRODUCT:	Cytomegalovirus (CMV) Monoclonal Antibody
IMMUNOGEN:	Affinity purified early antigen from MRC-5 infected cells infected with AD169 (ATCC).
FORM:	Ig purified from mouse Ascites in PBS, pH 7.6, 0.1% NaN <sub>3</sub> , sterile filtered.
SUBCLASS:	IgG <sub>2a</sub>
SPECIFICITY:	Immediate early nonstructural antigen (68-72kD). Does not cross react with EBV, ADV, VZV, or HSV.
STORAGE CONDITIONS:	4°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.