



PRODUCT SPECIFICATION

catalog number:	A28092
product description:	CMV pp65 Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	Specific for pp65 of CMV.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.