



PRODUCT SPECIFICATION

catalog number:	A29310
product description:	Calprotectin (S100A8/A9 Complex) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human calprotectin protein (native)
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
buffer:	15mM phosphate, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN ₃
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for human calprotectin S100A8/A9 complex.
pairing:	This antibody can be used as the detection antibody in a sandwich EIA and pairs with A29320 as the capture antibody.
notes/applications:	This antibody is known to work in latex agglutination and lateral flow assays.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.