



PRODUCT SPECIFICATION

catalog number:	A29350601P
product description:	Calprotectin Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Human Calprotectin Tetrameric Recombinant Protein
source:	Murine
form:	Liquid
purification method:	Protein A affinity chromatography
buffer:	PBS, pH 7.4
preservative:	0.09% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Recognizes the dimer S100A8 and S100A9
expiration date:	[lot specific]
storage:	-20°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.