



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A29351601P |
| product description: | Calprotectin Monoclonal Antibody |
| lot number: | [lot specific] |
| immunogen: | Human Calprotectin Tetrameric Recombinant Protein |
| source: | Murine |
| form: | Liquid |
| purification method: | Protein A affinity chromatography |
| buffer: | PBS, pH 7.4 |
| preservative: | 0.09% NaN ₃ |
| isotype: | IgG ₁ |
| concentration: | [lot specific] mg/mL |
| specificity: | Recognizes the dimer S100A8 and S100A9 |
| expiration date: | [lot specific] |
| storage: | -20°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |