



PRODUCT SPECIFICATION

CATALOG NUMBER: A32060

PRODUCT:	Gentamicin Monoclonal Antibody
PURIFICATION:	Protein A
PURITY:	>90% by SDS-PAGE
BUFFER:	10 mM Phosphate, 150 mM Na Cl, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁ , Kappa
STORAGE CONDITIONS:	2-8°C short term, -20 °C long term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.