



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A33140**

PRODUCT:	Human Ventricular Slow Skeletal Myosin Light Chain-1 (MLC-1) Monoclonal Antibody
IMMUNOGEN:	Purified MLC-1
PURIFICATION:	Protein G
BUFFER:	PBS, pH 7.3
PRESEVATIVE:	None
ISOTYPE:	IgG <sub>1kappa</sub>
STORAGE CONDITIONS:	-20°C
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.