



PRODUCT SPECIFICATION

catalog number:	A34012
product description:	HBsAg Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Highly purified HBsAg
source:	Murine ascites
form:	Liquid
purity:	>95% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN ₃
isotype:	IgG _{1, kappa}
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A34013, A34011 or G44020 (as the conjugate antibody).
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.