



PRODUCT SPECIFICATION

CATALOG NUMBER: A34060

PRODUCT:	HBsAg Monoclonal Antibody
FORM:	DEAE column chromatography purified immunoglobulin in 0.015 M KPO ₄ buffer, pH 7.40, containing 0.85% NaCl, 0.1% NaN ₃
AFFINITY CONSTANT:	2 x 10 ¹⁰ /mol
SUBCLASS:	IgG _{2b}
SPECIFICITY:	Reacts with the ad and ay subtype of the hepatitis B virus surface antigen.
STORAGE CONDITIONS:	2 - 8°C, Do not freeze
SUGGESTED USE:	For enzyme conjugation, radio labeling, or solid phase coating.
HAZARD/BIOHAZARD:	This material contains sodium azide a preservative, take appropriate precautions. There is no known biohazard associated with this product.