



PRODUCT SPECIFICATION

CATALOG NUMBER: A34090

PRODUCT:	HBsAg Monoclonal Antibody
FORM:	Protein G Sepharose purified in PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
SPECIFICITY:	Hepatitis B virus antigen immunoassay. Antibody reacts with the following HBsAg subtypes; ayw1, ayw2, ayw3, ayw4, ayr, adw2, adw4, adrq+, adrq-, ayw3 (Fer).
STORAGE CONDITIONS:	2 - 8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.