



PRODUCT SPECIFICATION

catalog number:	A34152
product description:	Hepatitis B Virus e Antigen (HBeAg) Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
purification method:	DEAE purified
purity:	A single band on Paragon electrophoresis gel
buffer:	0.01 M PBS, pH 7.4
preservative:	0.03% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	This antibody reacts with the hepatitis B virus e antigen
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A34162 (as the conjugate antibody).
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	2 - 8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.