



PRODUCT SPECIFICATION

catalog number:	A34171
product description:	Hepatitis B Virus e Antigen (HBeAg) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Hepatitis B virus
source:	Murine ascites
form:	Liquid
purity:	> 90%
buffer:	10mM PBS, 50% glycerin, pH 7.4
preservative:	0.09% NaN ₃
concentration:	[lot specific] mg/mL
specificity:	Reacts with hepatitis B virus e antigen.
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.