



PRODUCT SPECIFICATION

CATALOG NUMBER: A34213

PRODUCT:	Hepatitis B core antigen (HbcAg) Monoclonal Antibody
FORM:	DEAE purified Ig from mouse ascites
PURITY:	One single band on Paragon gel electrophoresis
BUFFER:	0.01 M PBS, pH 7.4
PRESERVATIVE:	0.03% sodium azide
SUBCLASS:	IgG ₃
STORAGE CONDITIONS:	2 - 8°C
SPECIFICITY:	Reacts with hepatitis B virus core antigen
PAIRING INFORMATION:	This mab pairs well with HBcAg mab A34214404P