



## PRODUCT SPECIFICATION

CATALOG NUMBER: A34215

PRODUCT DESCRIPTION:	Hepatitis B Virus core antigen (HbcAg) Monoclonal Antibody
FORM:	DEAE purified immunoglobulin from mouse ascites
PURITY:	One single band on Paragon gel electrophoresis
BUFFER:	0.01 M PBS, pH 7.4
PRESERVATIVE:	0.03% sodium azide
SPECIFICITY:	Reacts with hepatitis B virus core antigen
SUBCLASS:	IgG <sub>1</sub>
PAIRING:	This mab pairs well with HBcAg mab A34216 & A34217
STORAGE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.