



PRODUCT SPECIFICATION

CATALOG NUMBER: A34216

PRODUCT DESCRIPTION:	Hepatitis B Virus core antigen (HBcAg) Monoclonal Antibody
SPECIFICITY:	Reacts with hepatitis B virus core antigen
FORM:	DEAE purified immunoglobulin from mouse ascites
PURITY:	One single band on Paragon gel electrophoresis
BUFFER:	0.01 M PBS, pH 7.4
PRESERVATIVE:	0.03% sodium azide
SUBCLASS:	IgG _{2a}
PAIRING INFORMATION:	This MAb pairs well with HBcAg mab A34215
STORAGE:	2 - 8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.