



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A34220**

PRODUCT:	Hepatitis B virus core antigen (HBcAg) Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	Recombinant HBcAg.
FORM:	Each vial contains antibody in PBS, pH 7.4
PURIFICATION:	Purified by chromatography on protein G Sepharose. Purity is tested by electrophoresis.
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>2a</sub>
STORAGE CONDITIONS:	Store at 2-8°C.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.