



PRODUCT SPECIFICATION

CATALOG NUMBER: A34270

PRODUCT:	Glycogen Phosphorylase Isoenzyme BB (GPBB) Monoclonal Antibody
IMMUNOGEN:	Purified human GPBB.
PURITY:	Immunoglobulin purified by chromatography on protein A Sepharose. Purity is tested by electrophoresis.
BUFFER:	Antibody in PBS, pH 7.4
PRESERVATIVE:	0.1% of sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
PAIRING INFORMATION:	Mab #A34270 pairs with Affinity purified goat anti-GPBB peptide G #G-149-C.