



PRODUCT SPECIFICATION

CATALOG NUMBER: A34280

PRODUCT:	Glycogen Phosphorylase Isoenzyme BB (GPBB) Monoclonal Antibody
IMMUNOGEN:	Purified human GPBB.
PURIFICATION:	Immunoglobulin purified by chromatography on protein A sepharose.
BUFFER:	Antibody in PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG ₁
CROSS REACTIVITY:	No cross reaction with GPMM
STORAGE CONDITIONS:	2-8 ^o C