



PRODUCT SPECIFICATION

CATALOG NUMBER: A34290

PRODUCT:	Glycogen Phosphorylase Isoenzyme BB (GPBB) Monoclonal Antibody
SOURCE:	Immunoglobulin from mouse ascites
IMMUNOGEN:	Glycogen Phosphorylase Isoenzyme BB (GPBB)
FORM:	Purified antibody in PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
PURIFICATION:	Purified by chromatography on protein A Sepharose. Purity is tested electrophoresis.
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	4-8°C.
APPLICATION:	EIA, Western Blot