



PRODUCT SPECIFICATION

catalog number:	A34300
product description:	Hepatitis A Virus (HAV) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Cell culture
form:	Liquid
purification method:	Affinity chromatography
purity:	>70% pure by SDS-PAGE
buffer:	PBS
preservative:	0.01% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody recognizes the immunodominant site of HAV.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.