



PRODUCT SPECIFICATION

catalog number:	A34310
product description:	Hepatitis A Virus (HAV) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified HAV
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 6.8
preservative:	0.05% NaN ₃
isotype:	IgG _{2a} , kappa
concentration:	[lot specific] mg/mL
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.