



## PRODUCT SPECIFICATION

catalog number:	<b>A34501</b>
product description:	Cardiac Troponin I Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human cardiac troponin I subunit
form:	Liquid
purification method:	Protein A chromatography
buffer:	PBS, pH 7.4
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for aa 25-40 of human troponin I.
cross-reactivity:	There is no cross reactivity with human skeletal troponin I
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.