



PRODUCT SPECIFICATION

catalog number:	A34600
product description:	Troponin I Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Purified human cardiac troponin I subunit
form:	Liquid
purification method:	Protein A chromatography; 0.2µm filtered
buffer:	10mM sodium phosphate, 0.85% NaCl, pH 7.2
preservative:	0.05% NaN ₃
isotype:	IgG ₁ , kappa
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
manufacture date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.