



PRODUCT SPECIFICATION

CATALOG NUMBER: A34602

PRODUCT DESCRIPTION:	Troponin I (150-202) Monoclonal Antibody
IMMUNOGEN:	Troponin I
FORM:	Liquid
SOURCE:	Mouse ascites
PURITY:	>90% pure
BUFFER:	10 mM PBS, pH 7.4
PRESERVATIVE:	0.1% P300
APPLICATIONS:	Can be used in a sandwich ELISA as the detection antibody in combination with G-129-C (capture).
STORAGE:	2-8°C short term; -20°C long term
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.