



## PRODUCT SPECIFICATION

CATALOG NUMBER: A34603

PRODUCT DESCRIPTION:	Troponin I (162-202) Monoclonal Antibody
IMMUNOGEN:	Troponin I
FORM:	Liquid
SOURCE:	Mouse ascites
PURITY:	>90% pure
BUFFER:	10 mM PBS, pH 7.4
PRESERVATIVE:	0.1% P300
APPLICATIONS:	Can be used in a sandwich ELISA as in combination with A34600 (detection) and A34606 (coating).
STORAGE:	2-8°C short term; -20°C long term
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.