



PRODUCT SPECIFICATION

CATALOG NUMBER: A34604

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | Troponin I (190-196) Monoclonal Antibody |
| IMMUNOGEN: | Troponin I |
| FORM: | Liquid |
| SOURCE: | Mouse ascites |
| PURITY: | >90% pure |
| BUFFER: | 10 mM PBS, pH 7.4 |
| PRESERVATIVE: | 0.1% P300 |
| STORAGE: | 2-8°C short term; -20°C long term |
| HAZARD/BIOHAZARD: | There is no known biohazard associated with this product. |