



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A34630**

PRODUCT:	Troponin I Monoclonal Antibody
IMMUNOGEN:	Purified human cardiac Troponin I subunit.
FORM:	Liquid
PURIFICATION:	Immunoglobulin purified by protein A sepharose.
PURITY:	>95 % SDS-PAGE
BUFFER:	0.05M Tris/HCl, 0.1M NaCl, pH 7.2
PRESERVATIVE:	15 mM NaN <sub>3</sub>
SUBCLASS:	IgG <sub>2a</sub>
STORAGE CONDITIONS:	2-8 <sup>o</sup> C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.