



## PRODUCT SPECIFICATION

catalog number:	<b>A34708</b>
product description:	Troponin I Monoclonal Antibody, aa 185-196
lot number:	[lot specific]
source:	Murine; <i>in vitro</i>
form:	Liquid
purification method:	Protein A/G chromatography
purity:	≥95% by SDS-PAGE
buffer:	150 mM sodium chloride, 20 mM sodium phosphate, pH 7.4
preservative:	0.03% Proclin 300
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL (A280)
pairing:	This antibody can be used as the capture antibody, with A34707 as the detection.
expiration date:	[lot specific]
storage:	2-8°C short term (up to three months); -20°C long term (up to 36 months) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains proclin 300 as a preservative; take appropriate precaution. There is no known hazard/biohazard associated with this product.